

EVOLUTION PATIENT POSITIONERS

EVOLUTION Supine Positioners

REF	VSHP23	Supine Head Positioner Medium Style 2
REF	VSHP24	Supine Head Positioner Large Style 2
REF	VSKP13	Supine Knee Positioner Medium
REF	VSFP13	Supine Foot Positioner Medium

Intended for use on Operating Table Mattress to hold Patient stable during surgery.



Supine Foot Positioner



Supine Head Positioner

Supine Knee Positioner

VascoCare Medical has introduced innovative new Supine Positioners to the EVOLUTION Patient Positioner Range. These provide enhanced pressure relief with a proven anatomical design.

Features

- Unique and innovative design that **redistributes weight of the patient's body and provides enhanced pressure relief** compared to regular Patient Positioners.
- The positioners are made from two components. The inside of the product is comfortable **Visco-elastic foam**. The outer coating is a **seamless, non-porous Dual Coated Hi-Tech Membrane** with high tear strength making it a very robust product.
- The outer coating has **Excellent Abrasion Resistance** and is **Puncture/Burst Resistant**.
- It has excellent flexibility that **improves conformity to the skin**.
- Provides **Zero Interface Pressure** under the Bony Prominences.
- High strength bond between foam and coating creates an ultra strong laminated structure for a **long life product**.
- The outer coating provides an excellent **fluid-proof barrier without any seams** or other **bacteria harbouring features on the outer surface**.
- Completely **Latex free**.
- **Compatible with CT Scanners** as these do not cause artefacts.
- Head Positioners provide **excellent Neck Support and Head Stability** compared to regular Head Rings.
- Knee Positioners provide **excellent Support under the Leg** compared to regular Domed Patient Positioners.
- Foot Positioners provide **excellent Foot Support with Zero Pressure under the Heel** compared to regular Heel Pads.
- Competitively priced to make it a **cost effective option** against regular positioners.
- **Two Year Product Warranty**.



VascoCare Medical Limited
Enterprise Research Incubation Centre
Institute of Technology Carlow
Kilkenny Road, Carlow
Ireland
T: +353 (0)59 9175260
e: customercare@vascocare.com
w: www.vascocare.com

